Stable Micro Systems

Stress

g
Kg
N
mN

g/mm²
Kg/cm²
N/m²
Kg/m²
Pa
MPa

DR

KPa
Pa.s

°C

dB

g

Ohm
mOhm

formulas

Balance AED PT100
Define

See The Feeling

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

The diagram on the left shows the instrument's structure with three main components: 力臂 (Force Arm), 传动装置 (Transmission Device), and 底座 (Base). The photograph on the right shows the physical instrument, a Stable Micro Systems model, which is a white and blue benchtop device with a vertical force arm and a sample pan.

Lotun Science Best testing & Analysis Corp.

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The screenshot shows the software interface with a project list on the left and a graph of Force (g) vs. Time (s) on the right. The graph shows a single data point at approximately 10g. The interface includes a menu bar, a toolbar, and a status bar.

The top screenshot shows the '力量校正' (Force Calibration) step, where the user is prompted to place weights on the calibration platform and check the reading. The bottom screenshot shows the '高度校正' (Height Calibration) step, where the user is prompted to move the probe and check the reading. Both screenshots show the software interface with a graph and a list of parameters.

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Command	Program	Flags	Comment
Mark Value	1: Clear Graph Results		
Drop Anchor	2: Redraw		
Select Anchor	3: Search Forwards		
Search Forwards	4: Go to Minimum Time		
Search Backwards	5: Drop Anchor		
Go to Minimum Time	6: Set Threshold Force 5 g		
Go to Maximum Time	7: Go to Peak Top Force		
Set Threshold	8: F1 = Mark Value(Force (Curren... R		
Go to Peak Top	9: Drop Anchor		
Go to Peak Bottom	10: W1 = Area (Traditional)(Active...		
Area (Traditional)	11: Go to Force 5 g		
Gradient	12: Drop Anchor		
Mean	13: W2 = Area (Traditional)(Active...		
	14: R0 = W2 / W1	F	
	15: R = R0 * 100	F	
	16: Store Variable R as "彈性" "%"		
	17: Go to Absolute Maximum Valu...		
	18: F2 = Mark Value(Force (Curren... R		
	19: F = F2 - F1	F	
	20: Store Variable F as "韌性" "g"		
	21: Mark Value(Strain (%)) As 脆性 R		